



# Product Specification Summary

## PCR Amplification & Analysis

Catalog Number	40-5213-01
Product Name	PCR Master Mix (2X) 100 reactions (2 x 1.3 ml); 100 rxns
Size	100 rxns
Description	Ready to Use
Component/Note	
Component/Note	Add Primer and Template
Storage	Store at -20°C

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1056>



Product Manual [http://www.genelink.com/Literature/ps/M40-5213-XX\\_Ver3.2.pdf](http://www.genelink.com/Literature/ps/M40-5213-XX_Ver3.2.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	