



Product Specification Summary

PCR Amplification & Analysis

Catalog Number	40-5213-02
Product Name	PCR Master Mix (2X) 200 reactions (4 x 1.3 ml); 200 rxns
Size	200 rxns
Description	PCR Master Mix (2X); 200 reactions (4 x 1.3 ml)

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1057>



Product Manual http://www.genelink.com/Literature/ps/M40-5213-XX_Ver3.2.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	500 uL