



# Product Specification Summary

## PCR Amplification & Analysis

|                |  |
|----------------|--|
| Catalog Number | 40-5213-02   |
| Product Name   | PCR Master Mix (2X) 200 reactions (4 x 1.3 ml); 200 rxns |
| Size           | 200 rxns   |
| Description    | PCR Master Mix (2X); 200 reactions (4 x 1.3 ml)          |

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1057>



Product Manual [http://www.genelink.com/Literature/ps/M40-5213-XX\\_Ver3.2.pdf](http://www.genelink.com/Literature/ps/M40-5213-XX_Ver3.2.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

| Product                                  | Catalog No | Size   |
|--|------------|--------|
| Omni-Marker™ Universal unlabeled; 100 uL | 40-3005-01 | 100 uL |
| Omni-Marker™ Universal unlabeled; 1 mL   | 40-3005-10 | 1 mL   |
| Omni-Marker™ Low unlabeled; 500 uL       | 40-3006-05 | 500 uL |