



Product Specification Summary

GeneAssays™

Catalog Number	41-1036-01
Product Name	GeneAssays™ Genotyping PCR 2X Master Mix with ROX, 1 mL
Size	1 mL
Description	GeneAssays™ Genotyping PCR 2X Master Mix with ROX, 1 mL
Component/Note	GeneAssays™ Genotyping PCR 2X Master Mix with ROX, 1 mL

GeneAssays™ Genotyping PCR 2X Master Mix without dye and with passive ROX dye

The GeneAssays™ Genotyping PCR 2X Master Mix without dye and with passive ROX dye specifically formulated for qPCR, SNP and gene expression assays using fluorescent detection probes. The PCR master mix is validated for use with Gene Link designed and synthesized custom SNP and gene expression assays.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17348>



Product Manual http://www.genelink.com/Literature/ps/GeneAssays_2XMaster_Mix_V1.4.pdf



Product MSDS http://www.genelink.com/Literature/ps/MSDS_nonhaz_GeneAssays.pdf



Related Products

Product	Catalog No	Size
Custom SPCT (Synthetic Positive Control Template); up to 130mer; 500 pmols	26-6100-50	500 pmol
Custom SPCT (Synthetic Positive Control Template); 131-180mer; 500 pmols	26-6101-50	500 pmols

Custom SPCT (Synthetic Positive Control Template); 131-180mer; 2 nmols
190 Saw Mill River Road, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | www.genelink.com |
SS41-1036-01.pdf Print Date Version : May 28, 2017 Page 1 of 1

