



# Product Specification Summary

## GeneAssays™

Catalog Number 41-1121-02  
Product Name GeneAssays Custom SNP Genotyping Assay 20X Mix, Medium, 3,200X 10uL reactions  
Size 3200 Reactions  
Description GeneAssays Custom SNP Genotyping Assay 20X Mix, Medium, 3,200X 10uL reactions  
**GeneAssays Custom Gene SNP Genotyping Assay 20X Mix, Medium, 3,200X 10 uL reactions.**

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17351>



Product Manual [http://www.genelink.com/Literature/ps/GeneAssay\\_Custom\\_qPCR\\_Assay\\_PS\\_V3.2.pdf](http://www.genelink.com/Literature/ps/GeneAssay_Custom_qPCR_Assay_PS_V3.2.pdf)



Product MSDS [http://www.genelink.com/Literature/ps/MSDS\\_nonhaz\\_GeneAssays.pdf](http://www.genelink.com/Literature/ps/MSDS_nonhaz_GeneAssays.pdf)



### Related Products

Product	Catalog No	Size
Custom SPCT (Synthetic Positive Control Template); up to 130mer; 500 pmols	26-6100-50	500 pmol
Custom SPCT (Synthetic Positive Control Template); 131-180mer; 500 pmols	26-6101-50	500 pmols
Custom SPCT (Synthetic Positive Control Template); 131-180mer; 2 nmols	26-6101-02	2 nmols