



# Product Specification Summary

## GeneAssays™

Catalog Number 41-1141-02  
Product Name GeneAssays Custom Gene Expression Genotyping Assay 20X Mix, Medium, 3,200X 20 uL reactions  
Size 3200 Reactions  
Description GeneAssays Custom Gene Expression Genotyping Assay 20X Mix, Medium, 3,200X 20 uL reactions  
**GeneAssays Custom Gene Expression Genotyping Assay 20X Mix, Medium, 3,200X 20 uL reaction.**

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17354>



Product Manual [http://www.genelink.com/Literature/ps/GeneAssay\\_Custom\\_qPCR\\_Assay\\_PS\\_V3.2.pdf](http://www.genelink.com/Literature/ps/GeneAssay_Custom_qPCR_Assay_PS_V3.2.pdf)



Product MSDS [http://www.genelink.com/Literature/ps/MSDS\\_nonhaz\\_GeneAssays.pdf](http://www.genelink.com/Literature/ps/MSDS_nonhaz_GeneAssays.pdf)



### Related Products

Product	Catalog No	Size
Custom SPCT (Synthetic Positive Control Template); up to 130mer; 500 pmols	26-6100-50	500 pmol
Custom SPCT (Synthetic Positive Control Template); 131-180mer; 500 pmols	26-6101-50	500 pmols
Custom SPCT (Synthetic Positive Control Template); 131-180mer; 2 nmols	26-6101-02	2 nmols