



Product Specification Summary

GeneAssays™

Catalog Number 41-1141-03
Product Name GeneAssays Custom Gene Expression Genotyping Assay 20X Mix, Large, 6,400X 20 uL reactions
Size 6400 Reactions
Description GeneAssays Custom Gene Expression Genotyping Assay 20X Mix, Large, 6,400X 20 uL reactions
GeneAssays Custom Gene Expression Genotyping Assay 20X Mix, Large, 6,400X 20 uL reactions.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17355>



Product Manual http://www.genelink.com/Literature/ps/GeneAssay_Custom_qPCR_Assay_PS_V3.2.pdf



Product MSDS http://www.genelink.com/Literature/ps/MSDS_nonhaz_GeneAssays.pdf



Related Products

Product	Catalog No	Size
Custom SPCT (Synthetic Positive Control Template); up to 130mer; 500 pmols	26-6100-50	500 pmol
Custom SPCT (Synthetic Positive Control Template); 131-180mer; 500 pmols	26-6101-50	500 pmols
Custom SPCT (Synthetic Positive Control Template); 131-180mer; 2 nmols	26-6101-02	2 nmols