


Product Specification Summary

Storage Racks and Boxes

Catalog Number BB-3100-05-CL
Product Name Freezer Storage Box 100 Place for Microcentrifuge tubes; Clear, 5 each Pk
Size 1 Pk



Description Freezer Storage Box 100 Place for 2 mL micro tubes; 5X 10 Unit Pack

Component/Note Polypropylene 100 place box

Component/Note Freezer Storage Box 100 Place with hinged lid and alphanumeric numbering for 2 mL screw cap tubes and 1.5 mL screw cap tubes. Boxes made from virgin polypropylene withstands ultra low temperature freezers

Ultralow Freezer Storage Box. 100 Place

Durable freezer boxes that withstand ultra-low freezer temperatures of -80°C and manufactured using virgin polypropylene. Redesigned with sturdy lid with 3 hinges, alphanumeric numbered for convenient cataloging. Available in a selection of color for easy color coding.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17015>



Product Manual

http://www.genelink.com/Literature/ps/freezer-box_PS.pdf



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|---|---------------|--------|
| Microcentrifuge tube reversible rack Black; 5 each Pk | BB-3101-05-BK | 1 Pk |
| Microcentrifuge tube reversible rack Blue; 5 each Pk | BB-3101-05-BL | 1 Pk |
| Microcentrifuge tube reversible rack Assorted; 5 each Pk | BB-3101-05-AS | 1 Pk |
| Microcentrifuge tube reversible rack Green; 5 each Pk | BB-3101-05-GN | 1 Pk |
| Microcentrifuge tube reversible rack Natural; 5 each Pk | BB-3101-05-CL | 1 PK |
| Microcentrifuge tube reversible rack Red; 5 each Pk | BB-3101-05-RD | 1 Pk |
| Microcentrifuge tube reversible rack Yellow; 5 each Pk | BB-3101-05-YL | 1 Pk |
| Freezer Storage Box 100 Place for Microcentrifuge tubes; Blue, 1 each | BB-3100-01-BL | 1 Each |
| Freezer Storage Box 100 Place for Microcentrifuge tubes; Blue, 5 each Pk | BB-3100-05-BL | 1 Pk |
| Freezer Storage Box 100 Place for Microcentrifuge tubes; Blue, 2X 5 each Pk | BB-3100-10-BL | 2 Pk |
| Freezer Storage Box 100 Place for Microcentrifuge tubes; Assorted, 5 each Pk | BB-3100-05-AS | 1 Pk |
| Freezer Storage Box 100 Place for Microcentrifuge tubes; Yellow, 1 each | BB-3100-01-YL | 1 Each |
| Freezer Storage Box 100 Place for Microcentrifuge tubes; Yellow, 5 each Pk | BB-3100-05-YL | 1 Pk |
| Freezer Storage Box 100 Place for Microcentrifuge tubes; Yellow, 2X 5 each Pk | BB-3100-10-YL | 2 Pk |
| Freezer Storage Box 100 Place for Microcentrifuge tubes; Green, 1 each | BB-3100-01-GN | 1 Each |