


# Product Specification Summary

## BioMed Basics

# BioMed Basics

Catalog Number BB-3100-05-CL  
Product Name Freezer Storage Box 100 Place for Microcentrifuge tubes; Clear, 5 each Pk  
Size 1 Pk



Description Freezer Storage Box 100 Place for 2 mL micro tubes; 5X 10 Unit Pack

Component/Note Polypropylene 100 place box

Component/Note Freezer Storage Box 100 Place with hinged lid and alphanumeric numbering for 2 mL screw cap tubes and 1.5 mL screw cap tubes. Boxes made from virgin polypropylene withstands ultra low temperature freezers

### Ultralow Freezer Storage Box. 100 Place

Durable freezer boxes that withstand ultra-low freezer temperatures of  $-80^{\circ}\text{C}$  and manufactured using virgin polypropylene. Redesigned with sturdy lid with 3 hinges, alphanumeric numbered for convenient cataloging. Available in a selection of color for easy color coding.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17015>



Product Manual

[http://www.genelink.com/Literature/ps/freezer-box\\_PS.pdf](http://www.genelink.com/Literature/ps/freezer-box_PS.pdf)



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

| Product  | Catalog No    | Size     |
|--|---------------|----------|
| DNA Synthesis Column Luer Lock Set; Asstd. Colors      | GL-1010-50-ST | 50 set   |
| DNA Synthesis Column Standard Luer Lock; Natural Color | GL-1011-05-NC | 50 each  |
| DNA Synthesis Column Standard Luer Lock; Blue Color    | GL-1011-05-BL | 50 each  |
| DNA Synthesis Column Standard Luer Lock; Yellow Color  | GL-1011-05-YL | 50 each  |
| DNA Synthesis Column Insert with filter; Natural Color | GL-1012-05-NC | 50 each  |
| DNA Synthesis Column Screw Caps; Natural Color         | GL-1013-05-NC | 50 each  |
| DNA Synthesis Column Screw Caps; Green Color           | GL-1013-05-GN | 50 each  |
| DNA Synthesis Column Screw Caps; Red Color             | GL-1013-05-RD | 50 each  |
| DNA Synthesis Column Screw Caps; Blue Color            | GL-1013-05-BL | 50 each  |
| DNA Synthesis Column Screw Caps; Yellow Color          | GL-1013-05-YL | 50 each  |
| Micro tube clam shell holder 10 place; 100 Unit Pk     | GL-1002-10-NC | 1 Pk     |
| Micro tube clam shell holder 4 place; 10 Unit Pk       | GL-1001-01-NC | 1 Pk     |
| Micro tube clam shell holder 4 place; 50 Unit Pk       | GL-1001-05-NC | 1 Pk     |
| Micro tube clam shell holder 4 place; 100 Unit Pk      | GL-1001-10-NC | 1 Pk     |
| Cap Inserts for 2mL tube caps Green; 100 unit          | BB-4123-01-GN | 100 each |