

# Product Specification Summary

## Storage Racks and Boxes

Catalog Number BB-3101-05-BL  
Product Name Microcentrifuge tube reversible rack Blue; 5 each Pk  
Size 1 Pk



Description Microcentrifuge tube reversible rack for 0.5 mL, 1.5 mL and 2 mL tubes Blue  
Component/Note Microcentrifuge tube reversible rack for 0.5 mL, 1.5 mL and 2 mL tubes Blue  
Component/Note Microcentrifuge tube reversible rack for 0.5 mL, 1.5 mL and 2 mL tubes Blue

### **Microcentrifuge Tube Reversible Rack for 2mL , 1.5mL & 0.5mL Tubes**

Reversible rack for 0.5ml, 1.5ml and 2ml microcentrifuge tubes. Reversible for use with both 2ml and 1.5ml tubes on one side and 0.5ml tubes on the reverse side. A selection of colored racks is available

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17085>



Product Manual

<http://www.genelink.com/Literature/ps/none>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

| Product   | Catalog No    | Size   |
|---|---------------|--------|
| Microcentrifuge tube reversible rack Black; 5 each Pk                         | BB-3101-05-BK | 1 Pk   |
| Microcentrifuge tube reversible rack Blue; 5 each Pk                          | BB-3101-05-BL | 1 Pk   |
| Microcentrifuge tube reversible rack Assorted; 5 each Pk                      | BB-3101-05-AS | 1 Pk   |
| Microcentrifuge tube reversible rack Green; 5 each Pk                         | BB-3101-05-GN | 1 Pk   |
| Microcentrifuge tube reversible rack Natural; 5 each Pk                       | BB-3101-05-CL | 1 PK   |
| Microcentrifuge tube reversible rack Red; 5 each Pk                           | BB-3101-05-RD | 1 Pk   |
| Microcentrifuge tube reversible rack Yellow; 5 each Pk                        | BB-3101-05-YL | 1 Pk   |
| Freezer Storage Box 100 Place for Microcentrifuge tubes; Blue, 1 each         | BB-3100-01-BL | 1 Each |
| Freezer Storage Box 100 Place for Microcentrifuge tubes; Blue, 5 each Pk      | BB-3100-05-BL | 1 Pk   |
| Freezer Storage Box 100 Place for Microcentrifuge tubes; Blue, 2X 5 each Pk   | BB-3100-10-BL | 2 Pk   |
| Freezer Storage Box 100 Place for Microcentrifuge tubes; Assorted, 5 each Pk  | BB-3100-05-AS | 1 Pk   |
| Freezer Storage Box 100 Place for Microcentrifuge tubes; Yellow, 1 each       | BB-3100-01-YL | 1 Each |
| Freezer Storage Box 100 Place for Microcentrifuge tubes; Yellow, 5 each Pk    | BB-3100-05-YL | 1 Pk   |
| Freezer Storage Box 100 Place for Microcentrifuge tubes; Yellow, 2X 5 each Pk | BB-3100-10-YL | 2 Pk   |
| Freezer Storage Box 100 Place for Microcentrifuge tubes; Green, 1 each        | BB-3100-01-GN | 1 Each |