

# Product Specification Summary

## Storage Racks and Boxes

Catalog Number BB-3101-05-YL  
Product Name Microcentrifuge tube reversible rack Yellow; 5 each Pk  
Size 1 Pk



Description Microcentrifuge tube reversible rack for 0.5 mL, 1.5 mL and 2 mL tubes Yellow  
Component/Note Microcentrifuge tube reversible rack for 0.5 mL, 1.5 mL and 2 mL tubes Yellow  
Component/Note Microcentrifuge tube reversible rack for 0.5 mL, 1.5 mL and 2 mL tubes Yellow

### **Microcentrifuge Tube Reversible Rack for 2mL , 1.5mL & 0.5mL Tubes**

Reversible rack for 0.5ml, 1.5ml and 2ml microcentrifuge tubes. Reversible for use with both 2ml and 1.5ml tubes on one side and 0.5ml tubes on the reverse side. A selection of colored racks is available

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17087>



Product Manual

<http://www.genelink.com/Literature/ps/none>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
Microcentrifuge tube reversible rack Black; 5 each Pk	BB-3101-05-BK	1 Pk
Microcentrifuge tube reversible rack Blue; 5 each Pk	BB-3101-05-BL	1 Pk
Microcentrifuge tube reversible rack Assorted; 5 each Pk	BB-3101-05-AS	1 Pk
Microcentrifuge tube reversible rack Green; 5 each Pk	BB-3101-05-GN	1 Pk
Microcentrifuge tube reversible rack Natural; 5 each Pk	BB-3101-05-CL	1 PK
Microcentrifuge tube reversible rack Red; 5 each Pk	BB-3101-05-RD	1 Pk
Microcentrifuge tube reversible rack Yellow; 5 each Pk	BB-3101-05-YL	1 Pk
Freezer Storage Box 100 Place for Microcentrifuge tubes; Blue, 1 each	BB-3100-01-BL	1 Each
Freezer Storage Box 100 Place for Microcentrifuge tubes; Blue, 5 each Pk	BB-3100-05-BL	1 Pk
Freezer Storage Box 100 Place for Microcentrifuge tubes; Blue, 2X 5 each Pk	BB-3100-10-BL	2 Pk
Freezer Storage Box 100 Place for Microcentrifuge tubes; Assorted, 5 each Pk	BB-3100-05-AS	1 Pk
Freezer Storage Box 100 Place for Microcentrifuge tubes; Yellow, 1 each	BB-3100-01-YL	1 Each
Freezer Storage Box 100 Place for Microcentrifuge tubes; Yellow, 5 each Pk	BB-3100-05-YL	1 Pk
Freezer Storage Box 100 Place for Microcentrifuge tubes; Yellow, 2X 5 each Pk	BB-3100-10-YL	2 Pk
Freezer Storage Box 100 Place for Microcentrifuge tubes; Green, 1 each	BB-3100-01-GN	1 Each