

# Product Specification Summary

## 1.5 mL Microcentrifuge tubes QuickLock™ Cap

Catalog Number	BB-4110-05-AM
Product Name	Microcentrifuge tubes QuickLock™ Cap 1.5 mL Amber; 500 Each
Size	500 each



Description	Microcentrifuge tubes QuickLock™ Cap 1.5 mL Amber
Component/Note	Microcentrifuge tubes QuickLock™ Cap 1.5 mL Amber
Component/Note	Microcentrifuge tubes QuickLock™ Cap 1.5 mL Amber

### **1.5 mL QuickLock™ Snap Cap & 500 µL Standard Cap Microcentrifuge Tubes**

Microcentrifuge tubes are manufactured using virgin polypropylene and packaged under dust free conditions. All are RNase, DNase and pyrogen free. QuickLock™ cap design for microcentrifuge tubes are easy to open and close but the QuickLock™ mechanism prevents from accidental opening, withstands boiling aqueous solutions. Standard microcentrifuge tubes do not have the QuickLock™ snap cap mechanism. All microcentrifuge tubes are available in a selection of color for easy color coding.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17031>



Product Manual

<http://www.genelink.com/Literature/ps/non>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
Microcentrifuge tube reversible rack Black; 5 each Pk	BB-3101-05-BK	1 Pk
Microcentrifuge tube reversible rack Blue; 5 each Pk	BB-3101-05-BL	1 Pk
Microcentrifuge tube reversible rack Green; 5 each Pk	BB-3101-05-GN	1 Pk
Microcentrifuge tube reversible rack Natural; 5 each Pk	BB-3101-05-CL	1 PK
Microcentrifuge tube reversible rack Red; 5 each Pk	BB-3101-05-RD	1 Pk
Microcentrifuge tube reversible rack Yellow; 5 each Pk	BB-3101-05-YL	1 Pk
Microcentrifuge tubes QuickLock™ Cap 1.5 mL Amber; 500 Each	BB-4110-05-AM	500 each
Microcentrifuge tubes QuickLock™ Cap 1.5 mL Blue; 500 Each	BB-4110-05-BL	500 each
Microcentrifuge tubes QuickLock™ Cap 1.5 mL Gray; 500 Each	BB-4110-05-GY	500 each
Microcentrifuge tubes QuickLock™ Cap 1.5 mL Green; 500 Each	BB-4110-05-GN	500 each
Microcentrifuge tubes QuickLock™ Cap 1.5 mL Natural; 500 Each	BB-4110-05-NC	500 each
Microcentrifuge tubes QuickLock™ Cap 1.5 mL Orange; 500 Each	BB-4110-05-OR	500 each
Microcentrifuge tubes QuickLock™ Cap 1.5 mL Pink; 500 Each	BB-4110-05-PK	500 each
Microcentrifuge tubes QuickLock™ Cap 1.5 mL Violet; 500 Each	BB-4110-05-VT	500 Units
Microcentrifuge tubes QuickLock™ Cap 1.5 mL Yellow; 500 Each	BB-4110-05-YL	500 each