

# Product Specification Summary

## Microcentrifuge Tubes

Catalog Number	BB-4110-05-AM
Product Name	Microcentrifuge tubes QuickLock™ Cap 1.5 mL Amber; 500 Each
Size	500 each



Description	Microcentrifuge tubes QuickLock™ Cap 1.5 mL Amber
Component/Note	Microcentrifuge tubes QuickLock™ Cap 1.5 mL Amber
Component/Note	Microcentrifuge tubes QuickLock™ Cap 1.5 mL Amber

### **1.5 mL QuickLock™ Snap Cap & 500 µL Standard Cap Microcentrifuge Tubes**

Microcentrifuge tubes are manufactured using virgin polypropylene and packaged under dust free conditions. All are RNase, DNase and pyrogen free. QuickLock™ cap design for microcentrifuge tubes are easy to open and close but the QuickLock™ mechanism prevents from accidental opening, withstands boiling aqueous solutions. Standard microcentrifuge tubes do not have the QuickLock™ snap cap mechanism. All microcentrifuge tubes are available in a selection of color for easy color coding.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17031>



Product Manual

<http://www.genelink.com/Literature/ps/non>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
Cap Inserts for 2mL tube caps Green; 100 unit	BB-4123-01-GN	100 each
Cap Inserts for 2mL tube caps Red; 100 unit	BB-4123-01-RD	100 each
Cap Inserts for 2mL tube caps Violet; 100 unit	BB-4123-01-VT	100 each
Cap Inserts for 2mL tube caps White; 100 unit	BB-4123-01-WT	100 each
Cap Inserts for 2mL tube caps Yellow; 100 unit	BB-4123-01-YL	100 each
Cap Inserts for 2mL tube caps assorted colors; 100 unit	BB-4123-01-AS	100 each
Cap Inserts for 2mL tube caps Blue; 100 unit	BB-4123-01-BL	100 each
Cap Inserts for 2mL tube caps Orange; 100 unit	BB-4123-01-OR	100 each
Microcentrifuge conical screw cap tube 2mL Amber without caps; 500 Each	BB-4120-05-AM	500 each
Microcentrifuge conical screw cap tube 2mL Natural without caps; 500 Each	BB-4120-05-NC	500 each
Microcentrifuge conical tube 2mL with Amber screw cap; 500 Each	BB-4121-05-AM	500 each
Microcentrifuge conical tube 2mL with Blue screw cap; 500 Each	BB-4121-05-BL	500 each
Microcentrifuge conical tube 2mL with Gray screw cap; 500 Each	BB-4121-05-GY	500 each
Microcentrifuge conical tube 2mL with Green screw cap; 500 Each	BB-4121-05-GN	500 each
Microcentrifuge conical tube 2mL with Naturalscrew cap; 500 Each	BB-4121-05-NC	500 each