

Product Specification Summary

1.5 mL Microcentrifuge tubes QuickLock™ Cap

Catalog Number BB-4110-05-BL
Product Name Microcentrifuge tubes QuickLock™ Cap 1.5 mL Blue; 500 Each
Size 500 each



Description Microcentrifuge tubes 1.5 mL Blue
Component/Note Microcentrifuge tubes 1.5 mL Blue
Component/Note Microcentrifuge tubes 1.5 mL Blue

1.5 mL QuickLock™ Snap Cap & 500 µL Standard Cap Microcentrifuge Tubes

Microcentrifuge tubes are manufactured using virgin polypropylene and packaged under dust free conditions. All are RNase, DNase and pyrogen free. QuickLock™ cap design for microcentrifuge tubes are easy to open and close but the QuickLock™ mechanism prevents from accidental opening, withstands boiling aqueous solutions. Standard microcentrifuge tubes do not have the QuickLock™ snap cap mechanism. All microcentrifuge tubes are available in a selection of color for easy color coding.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17032>



Product Manual

<http://www.genelink.com/Literature/ps/non>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|---|---------------|-----------|
| Microcentrifuge tube reversible rack Black; 5 each Pk | BB-3101-05-BK | 1 Pk |
| Microcentrifuge tube reversible rack Blue; 5 each Pk | BB-3101-05-BL | 1 Pk |
| Microcentrifuge tube reversible rack Green; 5 each Pk | BB-3101-05-GN | 1 Pk |
| Microcentrifuge tube reversible rack Natural; 5 each Pk | BB-3101-05-CL | 1 PK |
| Microcentrifuge tube reversible rack Red; 5 each Pk | BB-3101-05-RD | 1 Pk |
| Microcentrifuge tube reversible rack Yellow; 5 each Pk | BB-3101-05-YL | 1 Pk |
| Microcentrifuge tubes QuickLock™ Cap 1.5 mL Amber; 500 Each | BB-4110-05-AM | 500 each |
| Microcentrifuge tubes QuickLock™ Cap 1.5 mL Blue; 500 Each | BB-4110-05-BL | 500 each |
| Microcentrifuge tubes QuickLock™ Cap 1.5 mL Gray; 500 Each | BB-4110-05-GY | 500 each |
| Microcentrifuge tubes QuickLock™ Cap 1.5 mL Green; 500 Each | BB-4110-05-GN | 500 each |
| Microcentrifuge tubes QuickLock™ Cap 1.5 mL Natural; 500 Each | BB-4110-05-NC | 500 each |
| Microcentrifuge tubes QuickLock™ Cap 1.5 mL Orange; 500 Each | BB-4110-05-OR | 500 each |
| Microcentrifuge tubes QuickLock™ Cap 1.5 mL Pink; 500 Each | BB-4110-05-PK | 500 each |
| Microcentrifuge tubes QuickLock™ Cap 1.5 mL Violet; 500 Each | BB-4110-05-VT | 500 Units |
| Microcentrifuge tubes QuickLock™ Cap 1.5 mL Yellow; 500 Each | BB-4110-05-YL | 500 each |