

Product Specification Summary

Microcentrifuge Tubes

Catalog Number BB-4123-01-WT
Product Name Cap Inserts for 2mL tube caps White; 100 unit
Size 100 each



Description Cap Inserts for 2mL tube caps; 100 unit White
Component/Note 100 unit White

2mL Screw Cap Conical & 2mL Screw Cap Skirted Self Standing Tubes

The 2ml screw cap microcentrifuge tubes are available as conical bottom without a skirt for easy viewing of contents at the bottom after centrifugation and with a skirt as a self-standing tube. Both are available in natural for regular use and amber color for light sensitive fluorescent contents. The screw caps are available with O-ring and another version with self-sealing cap; both with convenient attached tether to prevent capping mishaps. The caps are available in a selection of color for easy color coding.

Cap inserts for O-ring caps are available in 8 colors for cross color labeling and easy identification. These are manufactured using virgin polypropylene and packaged under dust free conditions. All are RNase, DNase and pyrogen free.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17023>



Product Manual

<http://www.genelink.com/Literature/ps/none>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|---|---------------|----------|
| Cap Inserts for 2mL tube caps Green; 100 unit | BB-4123-01-GN | 100 each |
| Cap Inserts for 2mL tube caps Red; 100 unit | BB-4123-01-RD | 100 each |
| Cap Inserts for 2mL tube caps Violet; 100 unit | BB-4123-01-VT | 100 each |
| Cap Inserts for 2mL tube caps White; 100 unit | BB-4123-01-WT | 100 each |
| Cap Inserts for 2mL tube caps Yellow; 100 unit | BB-4123-01-YL | 100 each |
| Cap Inserts for 2mL tube caps assorted colors; 100 unit | BB-4123-01-AS | 100 each |
| Cap Inserts for 2mL tube caps Blue; 100 unit | BB-4123-01-BL | 100 each |
| Cap Inserts for 2mL tube caps Orange; 100 unit | BB-4123-01-OR | 100 each |
| Microcentrifuge conical screw cap tube 2mL Amber without caps; 500 Each | BB-4120-05-AM | 500 each |
| Microcentrifuge conical screw cap tube 2mL Natural without caps; 500 Each | BB-4120-05-NC | 500 each |
| Microcentrifuge conical tube 2mL with Amber screw cap; 500 Each | BB-4121-05-AM | 500 each |
| Microcentrifuge conical tube 2mL with Blue screw cap; 500 Each | BB-4121-05-BL | 500 each |
| Microcentrifuge conical tube 2mL with Gray screw cap; 500 Each | BB-4121-05-GY | 500 each |
| Microcentrifuge conical tube 2mL with Green screw cap; 500 Each | BB-4121-05-GN | 500 each |
| Microcentrifuge conical tube 2mL with Naturalscrew cap; 500 Each | BB-4121-05-NC | 500 each |