

# Product Specification Summary

## Microcentrifuge Tubes

Catalog Number	BB-4130-05-AM
Product Name	Microcentrifuge selfstanding screw cap tube 2mL Amber
Size	500 Units



Description	Microcentrifuge self standing screw cap tube 2 mL Amber; tube only
Component/Note	Microcentrifuge self standing screw cap tube 2 mL Amber; tube only
Component/Note	Microcentrifuge self standing screw cap tube 2 mL Amber; tube only

### **2mL Screw Cap Conical & 2mL Screw Cap Skirted Self Standing Tubes**

The 2ml screw cap microcentrifuge tubes are available as conical bottom without a skirt for easy viewing of contents at the bottom after centrifugation and with a skirt as a self-standing tube. Both are available in natural for regular use and amber color for light sensitive fluorescent contents. The screw caps are available with O-ring and another version with self-sealing cap; both with convenient attached tether to prevent capping mishaps. The caps are available in a selection of color for easy color coding.

**Cap inserts for O-ring caps** are available in 8 colors for cross color labeling and easy identification. These are manufactured using virgin polypropylene and packaged under dust free conditions. All are RNase, DNase and pyrogen free.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17060>



Product Manual

<http://www.genelink.com/Literature/ps/none>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
Cap Inserts for 2mL tube caps Green; 100 unit	BB-4123-01-GN	100 each
Cap Inserts for 2mL tube caps Red; 100 unit	BB-4123-01-RD	100 each
Cap Inserts for 2mL tube caps Violet; 100 unit	BB-4123-01-VT	100 each
Cap Inserts for 2mL tube caps White; 100 unit	BB-4123-01-WT	100 each
Cap Inserts for 2mL tube caps Yellow; 100 unit	BB-4123-01-YL	100 each
Cap Inserts for 2mL tube caps assorted colors; 100 unit	BB-4123-01-AS	100 each
Cap Inserts for 2mL tube caps Blue; 100 unit	BB-4123-01-BL	100 each
Cap Inserts for 2mL tube caps Orange; 100 unit	BB-4123-01-OR	100 each
Microcentrifuge conical screw cap tube 2mL Amber without caps; 500 Each	BB-4120-05-AM	500 each
Microcentrifuge conical screw cap tube 2mL Natural without caps; 500 Each	BB-4120-05-NC	500 each
Microcentrifuge conical tube 2mL with Amber screw cap; 500 Each	BB-4121-05-AM	500 each
Microcentrifuge conical tube 2mL with Blue screw cap; 500 Each	BB-4121-05-BL	500 each
Microcentrifuge conical tube 2mL with Gray screw cap; 500 Each	BB-4121-05-GY	500 each
Microcentrifuge conical tube 2mL with Green screw cap; 500 Each	BB-4121-05-GN	500 each
Microcentrifuge conical tube 2mL with Naturalscrew cap; 500 Each	BB-4121-05-NC	500 each