

Product Specification Summary

PCR Tubes

Catalog Number BB-5130-05-GY
Product Name PCR Microtubes flat cap 0.2 mL Gray; 500 Each
Size 500 Each



Description PCR Microtubes flat cap 0.2 mL Gray; 500 Each
Component/Note PCR Microtubes flat cap 0.2 mL Gray; 500 Each
Component/Note PCR Microtubes flat cap 0.2 mL Gray; 500 Each

PCR QuickLock™ Tubes & PCR Standard Tubes

All PCR tubes are pre-sterilized and are RNase, DNase and pyrogen free. PCR tubes are manufactured using virgin polypropylene and packaged under dust free conditions. QuickLock™ cap design for PCR tubes are easy to open and close but the QuickLock™ mechanism prevents from accidental opening, withstands boiling aqueous solutions and is leak proof during thermal cycling. Standard PCR tubes do not have the QuickLock™ cap mechanism. All PCR tubes are available in a selection of color for easy color coding when setting up numerous reactions.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17207>



Product Manual

<http://www.genelink.com/Literature/ps/none>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|---|---------------|----------|
| PCR Microtubes flat cap 0.2 mL Amber; 500 Each | BB-5130-05-AM | 500 Each |
| PCR Microtubes flat cap 0.2 mL Blue; 500 Each | BB-5130-05-BL | 500 Each |
| PCR Microtubes flat cap 0.2 mL Gray; 500 Each | BB-5130-05-GY | 500 Each |
| PCR Microtubes flat cap 0.2 mL Green; 500 Each | BB-5130-05-GN | 500 Each |
| PCR Microtubes flat cap 0.2 mL Natural; 500 Each | BB-5130-05-NC | 500 Each |
| PCR Microtubes flat cap 0.2 mL Orange; 500 Each | BB-5130-05-OR | 500 Each |
| PCR Microtubes flat cap 0.2 mL Pink; 500 Each | BB-5130-05-PK | 500 Each |
| PCR Microtubes flat cap 0.2 mL Violet; 500 Each | BB-5130-05-VT | 500 Each |
| PCR Microtubes flat cap 0.2 mL Yellow; 500 Each | BB-5130-05-YL | 500 Each |
| PCR Microtubes QuickLock™ flat cap 0.2 mL Amber; 500 Each | BB-5140-05-AM | 500 Each |
| PCR Microtubes QuickLock™ flat cap 0.2 mL Blue; 500 Each | BB-5140-05-BL | 500 Each |
| PCR Microtubes QuickLock™ flat cap 0.2 mL Gray; 500 Each | BB-5140-05-GY | 500 Each |
| PCR Microtubes QuickLock™ flat cap 0.2 mL Green; 500 Each | BB-5140-05-GN | 500 Each |
| PCR Microtubes QuickLock™ flat cap 0.2 mL Natural; 500 Each | BB-5140-05-NC | 500 Each |
| PCR Microtubes QuickLock™ flat cap 0.2 mL Orange; 500 Each | BB-5140-05-OR | 500 Each |