

# Product Specification Summary

## PCR Tubes

Catalog Number BB-5130-05-PK  
Product Name PCR Microtubes flat cap 0.2 mL Pink; 500 Each  
Size 500 Each



Description PCR Microtubes flat cap 0.2 mL Pink; 500 Each  
Component/Note PCR Microtubes flat cap 0.2 mL Pink; 500 Each  
Component/Note PCR Microtubes flat cap 0.2 mL Pink; 500 Each

### PCR QuickLock™ Tubes & PCR Standard Tubes

All PCR tubes are pre-sterilized and are RNase, DNase and pyrogen free. PCR tubes are manufactured using virgin polypropylene and packaged under dust free conditions. QuickLock™ cap design for PCR tubes are easy to open and close but the QuickLock™ mechanism prevents from accidental opening, withstands boiling aqueous solutions and is leak proof during thermal cycling. Standard PCR tubes do not have the QuickLock™ cap mechanism. All PCR tubes are available in a selection of color for easy color coding when setting up numerous reactions.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17209>



Product Manual

<http://www.genelink.com/Literature/ps/none>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
PCR Microtubes flat cap 0.2 mL Amber; 500 Each	BB-5130-05-AM	500 Each
PCR Microtubes flat cap 0.2 mL Blue; 500 Each	BB-5130-05-BL	500 Each
PCR Microtubes flat cap 0.2 mL Gray; 500 Each	BB-5130-05-GY	500 Each
PCR Microtubes flat cap 0.2 mL Green; 500 Each	BB-5130-05-GN	500 Each
PCR Microtubes flat cap 0.2 mL Natural; 500 Each	BB-5130-05-NC	500 Each
PCR Microtubes flat cap 0.2 mL Orange; 500 Each	BB-5130-05-OR	500 Each
PCR Microtubes flat cap 0.2 mL Pink; 500 Each	BB-5130-05-PK	500 Each
PCR Microtubes flat cap 0.2 mL Violet; 500 Each	BB-5130-05-VT	500 Each
PCR Microtubes flat cap 0.2 mL Yellow; 500 Each	BB-5130-05-YL	500 Each
PCR Microtubes QuickLock™ flat cap 0.2 mL Amber; 500 Each	BB-5140-05-AM	500 Each
PCR Microtubes QuickLock™ flat cap 0.2 mL Blue; 500 Each	BB-5140-05-BL	500 Each
PCR Microtubes QuickLock™ flat cap 0.2 mL Gray; 500 Each	BB-5140-05-GY	500 Each
PCR Microtubes QuickLock™ flat cap 0.2 mL Green; 500 Each	BB-5140-05-GN	500 Each
PCR Microtubes QuickLock™ flat cap 0.2 mL Natural; 500 Each	BB-5140-05-NC	500 Each
PCR Microtubes QuickLock™ flat cap 0.2 mL Orange; 500 Each	BB-5140-05-OR	500 Each