

# Product Specification Summary

## BioMed Basics

# BioMed Basics

|                |                                                 |
|----------------|-------------------------------------------------|
| Catalog Number | BB-5130-05-YL                                   |
| Product Name   | PCR Microtubes flat cap 0.2 mL Yellow; 500 Each |
| Size           | 500 Each                                        |



|                |                                                 |
|----------------|-------------------------------------------------|
| Description    | PCR Microtubes flat cap 0.2 mL Yellow; 500 Each |
| Component/Note | PCR Microtubes flat cap 0.2 mL Yellow; 500 Each |
| Component/Note | PCR Microtubes flat cap 0.2 mL Yellow; 500 Each |

### PCR QuickLock™ Tubes & PCR Standard Tubes

All PCR tubes are pre-sterilized and are RNase, DNase and pyrogen free. PCR tubes are manufactured using virgin polypropylene and packaged under dust free conditions. QuickLock™ cap design for PCR tubes are easy to open and close but the QuickLock™ mechanism prevents from accidental opening, withstands boiling aqueous solutions and is leak proof during thermal cycling. Standard PCR tubes do not have the QuickLock™ cap mechanism. All PCR tubes are available in a selection of color for easy color coding when setting up numerous reactions.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17211>



Product Manual

<http://www.genelink.com/Literature/ps/none>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

| Product                                                | Catalog No    | Size     |
|--------------------------------------------------------|---------------|----------|
| DNA Synthesis Column Luer Lock Set; Asstd. Colors      | GL-1010-50-ST | 50 set   |
| DNA Synthesis Column Standard Luer Lock; Natural Color | GL-1011-05-NC | 50 each  |
| DNA Synthesis Column Standard Luer Lock; Blue Color    | GL-1011-05-BL | 50 each  |
| DNA Synthesis Column Standard Luer Lock; Yellow Color  | GL-1011-05-YL | 50 each  |
| DNA Synthesis Column Insert with filter; Natural Color | GL-1012-05-NC | 50 each  |
| DNA Synthesis Column Screw Caps; Natural Color         | GL-1013-05-NC | 50 each  |
| DNA Synthesis Column Screw Caps; Green Color           | GL-1013-05-GN | 50 each  |
| DNA Synthesis Column Screw Caps; Red Color             | GL-1013-05-RD | 50 each  |
| DNA Synthesis Column Screw Caps; Blue Color            | GL-1013-05-BL | 50 each  |
| DNA Synthesis Column Screw Caps; Yellow Color          | GL-1013-05-YL | 50 each  |
| Micro tube clam shell holder 10 place; 100 Unit Pk     | GL-1002-10-NC | 1 Pk     |
| Micro tube clam shell holder 4 place; 10 Unit Pk       | GL-1001-01-NC | 1 Pk     |
| Micro tube clam shell holder 4 place; 50 Unit Pk       | GL-1001-05-NC | 1 Pk     |
| Micro tube clam shell holder 4 place; 100 Unit Pk      | GL-1001-10-NC | 1 Pk     |
| Cap Inserts for 2mL tube caps Green; 100 unit          | BB-4123-01-GN | 100 each |