

Product Specification Summary

BioMed Basics

BioMed Basics

| | |
|----------------|--|
| Catalog Number | BB-5140-05-BL |
| Product Name | PCR Microtubes QuickLock™ flat cap 0.2 mL Blue; 500 Each |
| Size | 500 Each |



| | |
|----------------|--|
| Description | PCR Microtubes QuickLock™ flat cap 0.2 mL Blue; 500 Each |
| Component/Note | PCR Microtubes QuickLock™ flat cap 0.2 mL Blue; 500 Each |
| Component/Note | PCR Microtubes QuickLock™ flat cap 0.2 mL Blue; 500 Each |

PCR QuickLock™ Tubes & PCR Standard Tubes

All PCR tubes are pre-sterilized and are RNase, DNase and pyrogen free. PCR tubes are manufactured using virgin polypropylene and packaged under dust free conditions. QuickLock™ cap design for PCR tubes are easy to open and close but the QuickLock™ mechanism prevents from accidental opening, withstands boiling aqueous solutions and is leak proof during thermal cycling. Standard PCR tubes do not have the QuickLock™ cap mechanism. All PCR tubes are available in a selection of color for easy color coding when setting up numerous reactions.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17197>



Product Manual

<http://www.genelink.com/Literature/ps/none>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|--|---------------|----------|
| DNA Synthesis Column Luer Lock Set; Asstd. Colors | GL-1010-50-ST | 50 set |
| DNA Synthesis Column Standard Luer Lock; Natural Color | GL-1011-05-NC | 50 each |
| DNA Synthesis Column Standard Luer Lock; Blue Color | GL-1011-05-BL | 50 each |
| DNA Synthesis Column Standard Luer Lock; Yellow Color | GL-1011-05-YL | 50 each |
| DNA Synthesis Column Insert with filter; Natural Color | GL-1012-05-NC | 50 each |
| DNA Synthesis Column Screw Caps; Natural Color | GL-1013-05-NC | 50 each |
| DNA Synthesis Column Screw Caps; Green Color | GL-1013-05-GN | 50 each |
| DNA Synthesis Column Screw Caps; Red Color | GL-1013-05-RD | 50 each |
| DNA Synthesis Column Screw Caps; Blue Color | GL-1013-05-BL | 50 each |
| DNA Synthesis Column Screw Caps; Yellow Color | GL-1013-05-YL | 50 each |
| Micro tube clam shell holder 10 place; 100 Unit Pk | GL-1002-10-NC | 1 Pk |
| Micro tube clam shell holder 4 place; 10 Unit Pk | GL-1001-01-NC | 1 Pk |
| Micro tube clam shell holder 4 place; 50 Unit Pk | GL-1001-05-NC | 1 Pk |
| Micro tube clam shell holder 4 place; 100 Unit Pk | GL-1001-10-NC | 1 Pk |
| Cap Inserts for 2mL tube caps Green; 100 unit | BB-4123-01-GN | 100 each |