

## Product Specification Summary

### BioMed Basics

# BioMed Basics

Catalog Number	BB-5140-05-YL
Product Name	PCR Microtubes QuickLock™ flat cap 0.2 mL Yellow; 500 Each
Size	500 Each



Description	PCR Microtubes QuickLock™ flat cap 0.2 mL Yellow; 500 Each
Component/Note	PCR Microtubes QuickLock™ flat cap 0.2 mL Yellow; 500 Each
Component/Note	PCR Microtubes QuickLock™ flat cap 0.2 mL Yellow; 500 Each

### PCR QuickLock™ Tubes & PCR Standard Tubes

All PCR tubes are pre-sterilized and are RNase, DNase and pyrogen free. PCR tubes are manufactured using virgin polypropylene and packaged under dust free conditions. QuickLock™ cap design for PCR tubes are easy to open and close but the QuickLock™ mechanism prevents from accidental opening, withstands boiling aqueous solutions and is leak proof during thermal cycling. Standard PCR tubes do not have the QuickLock™ cap mechanism. All PCR tubes are available in a selection of color for easy color coding when setting up numerous reactions.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17212>



Product Manual

<http://www.genelink.com/Literature/ps/NONE>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
DNA Synthesis Column Luer Lock Set; Asstd. Colors	GL-1010-50-ST	50 set
DNA Synthesis Column Standard Luer Lock; Natural Color	GL-1011-05-NC	50 each
DNA Synthesis Column Standard Luer Lock; Blue Color	GL-1011-05-BL	50 each
DNA Synthesis Column Standard Luer Lock; Yellow Color	GL-1011-05-YL	50 each
DNA Synthesis Column Insert with filter; Natural Color	GL-1012-05-NC	50 each
DNA Synthesis Column Screw Caps; Natural Color	GL-1013-05-NC	50 each
DNA Synthesis Column Screw Caps; Green Color	GL-1013-05-GN	50 each
DNA Synthesis Column Screw Caps; Red Color	GL-1013-05-RD	50 each
DNA Synthesis Column Screw Caps; Blue Color	GL-1013-05-BL	50 each
DNA Synthesis Column Screw Caps; Yellow Color	GL-1013-05-YL	50 each
Micro tube clam shell holder 10 place; 100 Unit Pk	GL-1002-10-NC	1 Pk
Micro tube clam shell holder 4 place; 10 Unit Pk	GL-1001-01-NC	1 Pk
Micro tube clam shell holder 4 place; 50 Unit Pk	GL-1001-05-NC	1 Pk
Micro tube clam shell holder 4 place; 100 Unit Pk	GL-1001-10-NC	1 Pk
Cap Inserts for 2mL tube caps Green; 100 unit	BB-4123-01-GN	100 each