

Product Specification Summary

Fluorescent Molecular Probes

Catalog Number	FMPPM
Product Name	Fluorescent Molecular Probes - Product Manual
Size	1 Doc
Description	Fluorescent Molecular Probes - Product Manual
Component/Note	Fluorescent Molecular Probes - Product Manual
Component/Note	Fluorescent Molecular Probes - Product Manual
Fluorescent Molecular Probes - Product Manual	

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=948>



Product Manual <http://www.genelink.com/Literature/ps/PG-MB-Ver5.3.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
SmartBase™ Fluorescent Molecular Probes	SB-FMP	1 each
QPCR_Manual	M_QPCR	1 each
5'-Fam/3'-Quencher 10 umol scale; 600 nmol yield	26-8121-10	600 nmols