

Product Specification Summary

Clam Shell Tube Box

Catalog Number	GL-1001-05-NC
Product Name	Micro tube clam shell holder 4 place; 50 Unit Pk
Size	1 Pk



Description	4 Place PVC Clam Shell
Component/Note	Clear with no pigmentation

Clamshell Microtube Shipping Container

PVC Clamshell Holders 4 Place & 10 Place for Storage and Shipping. Custom designed specifically for 2mL microcentrifuge tube shipping and storage. The tube slots securely clamp the tube and will not dislodge or fall off even when opened upside down. Easy clamshell lid closes firmly. Containers can also be used for convenient freezer storage.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17017>



Product Manual

http://www.genelink.com/Literature/ps/freezer-box_PS.pdf



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Micro tube clam shell holder 10 place; 100 Unit Pk	GL-1002-10-NC	1 Pk
Micro tube clam shell holder 4 place; 10 Unit Pk	GL-1001-01-NC	1 Pk
Micro tube clam shell holder 4 place; 50 Unit Pk	GL-1001-05-NC	1 Pk
Micro tube clam shell holder 4 place; 100 Unit Pk	GL-1001-10-NC	1 Pk
Micro tube clam shell holder 10 place; 50 Unit Pk	GL-1002-05-NC	1 Pk
Micro tube clam shell holder 10 place; 10 Unit Pk	GL-1002-01-NC	1 Pk