



Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17020>



Product Manual

[http://www.genelink.com/Literature/ps/freezer-box\\_PS.pdf](http://www.genelink.com/Literature/ps/freezer-box_PS.pdf)



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
DNA Synthesis Column Luer Lock Set; Asstd. Colors	GL-1010-50-ST	50 set
DNA Synthesis Column Standard Luer Lock; Natural Color	GL-1011-05-NC	50 each
DNA Synthesis Column Standard Luer Lock; Blue Color	GL-1011-05-BL	50 each
DNA Synthesis Column Standard Luer Lock; Yellow Color	GL-1011-05-YL	50 each
DNA Synthesis Column Insert with filter; Natural Color	GL-1012-05-NC	50 each
DNA Synthesis Column Screw Caps; Natural Color	GL-1013-05-NC	50 each
DNA Synthesis Column Screw Caps; Green Color	GL-1013-05-GN	50 each
DNA Synthesis Column Screw Caps; Red Color	GL-1013-05-RD	50 each
DNA Synthesis Column Screw Caps; Blue Color	GL-1013-05-BL	50 each
DNA Synthesis Column Screw Caps; Yellow Color	GL-1013-05-YL	50 each
Micro tube clam shell holder 10 place; 100 Unit Pk	GL-1002-10-NC	1 Pk
Micro tube clam shell holder 4 place; 10 Unit Pk	GL-1001-01-NC	1 Pk
Micro tube clam shell holder 4 place; 50 Unit Pk	GL-1001-05-NC	1 Pk
Micro tube clam shell holder 4 place; 100 Unit Pk	GL-1001-10-NC	1 Pk
Cap Inserts for 2mL tube caps Green; 100 unit	BB-4123-01-GN	100 each