

Product Specification Summary

Fluorescent Molecular Probes

Catalog Number M_QPCR
Product Name QPCR_Manual
Size 1 each
Description QPCR_Manual

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17180>



Product Manual http://www.genelink.com/Literature/ps/Control_RT_Probes_Manual_V2.3.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
SmartBase™ Fluorescent Molecular Probes	SB-FMP	1 each
QPCR_Manual	M_QPCR	1 each
5'-Fam/3'-Quencher 10 umol scale; 600 nmol yield	26-8121-10	600 nmols