



Product Specifications


Custom Oligo Synthesis, antisense oligos, RNA oligos, chimeric oligos, Fluorescent dyes, Affinity Ligands, Spacers & Linkers, Duplex Stabilizers, Minor bases, labeled oligos, Molecular Beacons, siRNA, phosphonates Locked Nucleic Acids (LNA); 2'-5' linked Oligos

Oligo Modifications

For research use only. Not for use in diagnostic procedures for clinical purposes.

vP-mU

Category	Antisense & siRNA
Modification Code	vP-mU
Reference Catalog Number	26-6597
5 Prime	Y
3 Prime	N
Internal	N
Molecular Weight(mw)	396.18



vinyl phosphonate 2'O methyl U [vP-mU]
[26-6497-XX]

vP-mU modification has a setup charge of \$350.00 per order for special synthesis reagents .

Antisense Oligos (ODN) & siRNA Oligo Modifications

Click here for more information on antisense modifications, design & applications.