


Product Specification

2 mL Tube Freezer Storage Boxes, 2 mL Tube Shipping Containers
 Reversible 2 mL & 0.5 mL racks, Screw cap tubes
 1.5 mL and 0.5 mL micro centrifuge tubes
 Pipet tips and racks, PCR tubes

Gene Link introduces quality plastic products in partnership with Biomed Basics



Freezer Storage Box

Freezer Storage Box 100 Place with hinged lid and alphanumeric numbering for 2 mL screw cap tubes and 1.5 mL snap cap tubes. Boxes made from virgin polypropylene withstands ultra low temperature freezers.

	Catalog No.	Product Description	Unit Size	Price
	BB-3100-01-CL	Freezer Storage Box 100 Place for 2 mL micro tubes; 1 Unit Pack	1	\$6.25
	BB-3100-05-CL	Freezer Storage Box 100 Place for 2 mL micro tubes; 5 Unit Pack	5	\$27.60
	BB-3100-10-CL	Freezer Storage Box 100 Place for 2 mL micro tubes; 10 Unit Pack	10	\$49.50

Micro Tube Shipping Container

Clam shell holder PVC 4 place and 10 place containers for storage and secure transportation.

	Catalog No.	Product Description	Unit Size	Price
	GL-1001-01-NC	Micro tube clam shell holder 4 place; 10 Unit Pack	10	\$3.75
	GL-1001-05-NC	Micro tube clam shell holder 4 place; 50 Unit Pack	50	\$16.50
	GL-1001-10-NC	Micro tube clam shell holder 4 place; 100 Unit Pack	100	\$29.75
	GL-1002-01-NC	Micro tube clam shell holder 10 place; 10 Unit Pack	10	\$6.25
	GL-1002-05-NC	Micro tube clam shell holder 10 place; 50 Unit Pack	50	\$27.00
	GL-1002-10-NC	Micro tube clam shell holder 10 place; 100 Unit Pack	100	\$48.50

BioMed Basics

 **Gene Link™**