



Product Specification Summary

Myotonic Dystrophy

Catalog Number	40-2026-11
Product Name	Myotonic Dystrophy Genemer Kit; 1 Kit
Size	1 Kit
Description	Myotonic Dystrophy Genotyping
Component/Note	CTG Repeat Amplification

Genemer™ are optimized primer pairs for amplification of the gene fragment of interest, particularly those spanning a mutation. A wide range is available from single base pair mutation sites, the complex triple repeat disorders and various pathogens. Genemer™ kits are complete, easy-to-use systems for reliable genotyping of large repeats in certain triple repeat disorder. The Genemer™ product contains 10 nmols of the primer pair that is sufficient for 400 regular 50 µl amplifications. The Genemer™ Kit contains reagents as well as control DNA sufficient for 100 amplifications.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=1024>



Product Manual

http://www.genelink.com/Literature/ps/M40-2026-11_Ver5.4.pdf



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Myotonic Dystrophy GeneProber GLDM2 (Pst I) unlabeled Probe; 500 ng	40-2026-39	500 ng
Myotonic Dystrophy GScan™ TET Kit	40-2026-15TT	1 Kit
Myotonic Dystrophy GScan™ HEX Kit	40-2026-15HX	1 Kit
Myotonic Dystrophy GScan™ Cy3 Kit	40-2026-15C3	1 Kit
Myotonic Dystrophy GScan™ Cy5 Kit	40-2026-15C5	1 Kit
Myotonic Dystrophy GeneProber GLDM3 unlabeled Probe; 500 ng	40-2026-38	500 ng
Myotonic Dystrophy GScan™ 6-FAM Kit; 1 Kit	40-2026-15FM	1 Kit
Myotonic Dystrophy GScan™ 6-FAM Kit (20rxns); 1 Kit	40-2026-15FMS	1 Kit
Myotonic Dystrophy Genemer Kit; 1 Kit	40-2026-11	1 Kit
DM CTG-6-MB-CY3; 25 ug	45-0210-06	25 ug
DM CTG-6-MB-CY5; 25 ug	45-0220-06	25 ug
DM CTG-6-MB-FAM; 25 ug	45-0230-06	25 ug
DM CTG-6-MB-FI; 25 ug	45-0240-06	25 ug
DM CTG-6-MB-HEX; 25 ug	45-0250-06	25 ug
DM CTG-6-MB-TET; 25 ug	45-0260-06	25 ug