



# Product Specification Summary

## PCR Amplification & Analysis

Catalog Number 40-5200-30  
Product Name Taq DNA Polymerase 300 units 5 u/ul 60 uL; 300 units  
Size 300 units  
Description Taq DNA Polymerase 300 units 5 u/ul 60 uL; 300 units  
Component/Note 5 u/ul  
Taq DNA Polymerase 300 units 5 u/ul 60 uL; 300 units

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1058>



Product Manual [http://www.genelink.com/Literature/ps/M40-5200-XX\\_Ver3.1.pdf](http://www.genelink.com/Literature/ps/M40-5200-XX_Ver3.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	500 uL