



Product Specification Summary

Sequencing Primers

Catalog Number 26-3000-12
Product Name SP6 Universal 18mer; 25 ug
Size 25 ug

This SP6 RNA polymerase primer is complementary to the conserved class class III sequence. Sequencing data will be obtained for all templates containing the Class III T7 RNA polymerase promoter sequences. The product is supplied as a lyophilized powder, after reconstitution store at -20oC. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=153>



Product Manual http://www.genelink.com/Literature/ps/T7_T3_seq_primer_PSV3.1.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|--------------------|------------|-------|
| M13 Rev Seq; 25 ug | 26-3000-52 | 25 ug |
| M13 For Seq; 25 ug | 26-3000-51 | 25 ug |
| T7 Seq Long | 26-3000-55 | 25 ug |