



# Product Specification Summary

## Genetic Tools & Reagents

Catalog Number	40-3015-50
Product Name	Agarose HiRes Ultra Pure Molecular Biology Grade; 500 g
Size	500 g
Description	Agarose HiRes Ultra Pure Molecular Biology Grade; 500 gms
Component/Note	Agarose HiRes is certified Ultra Pure molecular biology grade DNase and RNase-free agarose powder. It is specifically recommended for resolution of short fragments ranging in size between 20 and 800 bp
Component/Note	Agarose HiRes is certified Ultra Pure molecular biology grade DNase and RNase-free agarose powder. It is specifically recommended for resolution of short fragments ranging in size between 20 bp and 800 bp and is an excellent substitute for polyacrylamide electrophoresis for resolution of short DNA fragments. HiRes Agarose is commonly used for electrophoretic resolution of fragments obtained from amplification of short tandem repeats (STR's), di, tri and tetra-nucleotide repeats, and other polymorphic loci.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=16988>



Product Manual

<http://www.genelink.com/Literature/ps/PS40-3015-XXV1.2.pdf>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
Omni-Array™ mRN Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni- Marker™ Universal unlabeled; 500 uL	40-3005-05	500 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni- Marker™ Low unlabeled; 100 uL	40-3006-01	100 uL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	500 uL
Omni- Marker™ Low unlabeled; 1 mL	40-3006-10	1 mL
GScan Marker Tamra Labeled 50 bp - 600 bp; 100 uL	40-3061-01	100 uL
GScan Marker Tamra Labeled 50 bp - 600 bp; 500 uL	40-3061-05	500 uL
Omni-Marker™ GScan-1 Tet labeled 50 bp-1kb	40-3071-01	100 ul
Omni-Marker™ GScan-1 Tet labeled 50 bp-1kb	40-3071-05	500 ul
Omni-Marker™ GScan-2 Tet labeled 50 bp- 600 bp	40-3072-01	100 ul
Omni-Marker™ GScan-2 Tet labeled 50 bp- 600 bp	40-3072-05	500 ul
GScan Marker Hex Labeled 50 bp - 600 bp; 100 uL	40-3081-01	100 uL