

Product Specification Summary

BioMed Basics

BioMed Basics

| | |
|----------------|---|
| Catalog Number | BB-4122-05-OR |
| Product Name | Screw Caps Orange with O-ring Tethered for 2 mL microcentrifuge tubes; 500 Each |
| Size | 500 each |



| | |
|----------------|---|
| Description | Screw Caps Tethered for 2 mL microcentrifuge tubes with O-ring Orange |
| Component/Note | Screw Caps Tethered for 2 mL microcentrifuge tubes with O-ring Orange |
| Component/Note | Screw Caps Tethered for 2 mL microcentrifuge tubes with O-ring Orange |

2mL Screw Cap Conical & 2mL Screw Cap Skirted Self Standing Tubes

The 2ml screw cap microcentrifuge tubes are available as conical bottom without a skirt for easy viewing of contents at the bottom after centrifugation and with a skirt as a self-standing tube. Both are available in natural for regular use and amber color for light sensitive fluorescent contents. The screw caps are available with O-ring and another version with self-sealing cap; both with convenient attached tether to prevent capping mishaps. The caps are available in a selection of color for easy color coding.

Cap inserts for O-ring caps are available in 8 colors for cross color labeling and easy identification. These are manufactured using virgin polypropylene and packaged under dust free conditions. All are RNase, DNase and pyrogen free.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17077>



Product Manual

<http://www.genelink.com/Literature/ps/none>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|--|---------------|----------|
| DNA Synthesis Column Luer Lock Set; Asstd. Colors | GL-1010-50-ST | 50 set |
| DNA Synthesis Column Standard Luer Lock; Natural Color | GL-1011-05-NC | 50 each |
| DNA Synthesis Column Standard Luer Lock; Blue Color | GL-1011-05-BL | 50 each |
| DNA Synthesis Column Standard Luer Lock; Yellow Color | GL-1011-05-YL | 50 each |
| DNA Synthesis Column Insert with filter; Natual Color | GL-1012-05-NC | 50 each |
| DNA Synthesis Column Screw Caps; Natural Color | GL-1013-05-NC | 50 each |
| DNA Synthesis Column Screw Caps; Green Color | GL-1013-05-GN | 50 each |
| DNA Synthesis Column Screw Caps; Red Color | GL-1013-05-RD | 50 each |
| DNA Synthesis Column Screw Caps; Blue Color | GL-1013-05-BL | 50 each |
| DNA Synthesis Column Screw Caps; Yellow Color | GL-1013-05-YL | 50 each |
| Micro tube clam shell holder 10 place; 100 Unit Pk | GL-1002-10-NC | 1 Pk |
| Micro tube clam shell holder 4 place; 10 Unit Pk | GL-1001-01-NC | 1 Pk |
| Micro tube clam shell holder 4 place; 50 Unit Pk | GL-1001-05-NC | 1 Pk |
| Micro tube clam shell holder 4 place; 100 Unit Pk | GL-1001-10-NC | 1 Pk |
| Cap Inserts for 2mL tube caps Green; 100 unit | BB-4123-01-GN | 100 each |