

Product Specification Summary

BioMed Basics

BioMed Basics

Catalog Number	BB-4122-05-YL
Product Name	Screw Caps Yellow with O-ring Tethered for 2 mL microcentrifuge tubes; 500 Each
Size	500 each



Description	Screw Caps Tethered for 2 mL microcentrifuge tubes with O-ring Yellow
Component/Note	Screw Caps Tethered for 2 mL microcentrifuge tubes with O-ring Yellow
Component/Note	Screw Caps Tethered for 2 mL microcentrifuge tubes with O-ring Yellow

2mL Screw Cap Conical & 2mL Screw Cap Skirted Self Standing Tubes

The 2ml screw cap microcentrifuge tubes are available as conical bottom without a skirt for easy viewing of contents at the bottom after centrifugation and with a skirt as a self-standing tube. Both are available in natural for regular use and amber color for light sensitive fluorescent contents. The screw caps are available with O-ring and another version with self-sealing cap; both with convenient attached tether to prevent capping mishaps. The caps are available in a selection of color for easy color coding.

Cap inserts for O-ring caps are available in 8 colors for cross color labeling and easy identification. These are manufactured using virgin polypropylene and packaged under dust free conditions. All are RNase, DNase and pyrogen free.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17080>



Product Manual

<http://www.genelink.com/Literature/ps/none>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
DNA Synthesis Column Luer Lock Set; Asstd. Colors	GL-1010-50-ST	50 set
DNA Synthesis Column Standard Luer Lock; Natural Color	GL-1011-05-NC	50 each
DNA Synthesis Column Standard Luer Lock; Blue Color	GL-1011-05-BL	50 each
DNA Synthesis Column Standard Luer Lock; Yellow Color	GL-1011-05-YL	50 each
DNA Synthesis Column Insert with filter; Natual Color	GL-1012-05-NC	50 each
DNA Synthesis Column Screw Caps; Natural Color	GL-1013-05-NC	50 each
DNA Synthesis Column Screw Caps; Green Color	GL-1013-05-GN	50 each
DNA Synthesis Column Screw Caps; Red Color	GL-1013-05-RD	50 each
DNA Synthesis Column Screw Caps; Blue Color	GL-1013-05-BL	50 each
DNA Synthesis Column Screw Caps; Yellow Color	GL-1013-05-YL	50 each
Micro tube clam shell holder 10 place; 100 Unit Pk	GL-1002-10-NC	1 Pk
Micro tube clam shell holder 4 place; 10 Unit Pk	GL-1001-01-NC	1 Pk
Micro tube clam shell holder 4 place; 50 Unit Pk	GL-1001-05-NC	1 Pk
Micro tube clam shell holder 4 place; 100 Unit Pk	GL-1001-10-NC	1 Pk
Cap Inserts for 2mL tube caps Green; 100 unit	BB-4123-01-GN	100 each