

# Product Specification Summary

## Fluorescent Molecular Probes

Catalog Number M\_QPCR  
Product Name QPCR\_Manual  
Size 1 each  
Description QPCR\_Manual

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17180>



Product Manual [http://www.genelink.com/Literature/ps/Control\\_RT\\_Probes\\_PG\\_V3.2.pdf](http://www.genelink.com/Literature/ps/Control_RT_Probes_PG_V3.2.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
SmartBase™ Fluorescent Molecular Probes	SB-FMP	1 each
QPCR_Manual	M_QPCR	1 each
Fluorescent Probe, 5'-Fam/3'-Quencher, 600 nmol purified yield	26-8121-10	600 nmols