



# Product Specification Summary

## PCR Amplification & Analysis

|                |                                 |
|----------------|---------------------------------|
| Catalog Number | 40-3000-16                      |
| Product Name   | DEPC Treated Water; 10 X 1.6 mL |
| Size           | 1.6 mL                          |
| Description    | DEPC Treated Water; 10 X 1.6 mL |

Gene Link DEPC treated water and Nuclease Free water (DEPC free) has been treated and autoclaved water that are certified to be free of nucleases and nucleic acids. Nuclease Free water (DEPC free) has not been treated with DEPC and is certified for all molecular biology applications requiring nuclease and nucleic acid free water. The smaller 1.6 mL product size is specifically recommended for RT-PCR (QPCR) applications. Quality Control includes RT-PCR, standard PCR and RNA gels without any degradation certifies the product to free of non-specific endonuclease, exonuclease and RNase activity.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17181>



Product Manual [http://www.genelink.com/Literature/ps/PS40-30-DNA&RNA\\_Soln\\_Ver2.1.pdf](http://www.genelink.com/Literature/ps/PS40-30-DNA&RNA_Soln_Ver2.1.pdf)



Product MSDS [http://www.genelink.com/Literature/ps/SDS\\_40-DNA&RNA\\_solutions.pdf](http://www.genelink.com/Literature/ps/SDS_40-DNA&RNA_solutions.pdf)



### Related Products

| Product                                  | Catalog No | Size   |
|--|------------|--------|
| Omni-Marker™ Universal unlabeled; 100 uL | 40-3005-01 | 100 uL |
| Omni-Marker™ Universal unlabeled; 1 mL   | 40-3005-10 | 1 mL   |
| Omni-Marker™ Low unlabeled; 500 uL       | 40-3006-05 |        |