



Product Specification Summary

GeneAssays™

Catalog Number 41-1121-03
Product Name GeneAssays Custom SNP Genotyping Assay 20X Mix, Large, 6,400X 10uL reactions
Size 6400 Reactions
Description GeneAssays Custom SNP Genotyping Assay 20X Mix, Large, 6,400X 10µL reactions
GeneAssays™ Custom SNP Genotyping Assay 20X Mix, Large, 6,400X 10 uL reactions.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17352>



Product Manual http://www.genelink.com/Literature/ps/GeneAssay_Custom_qPCR_Assay_PS_V3.2.pdf



Product MSDS http://www.genelink.com/Literature/ps/MSDS_nonhaz_GeneAssays.pdf



Related Products

Product	Catalog No	Size
Custom SPCT (Synthetic Positive Control Template); up to 130mer; 500 pmols	26-6100-50	500 pmol
Custom SPCT (Synthetic Positive Control Template); 131-180mer; 500 pmols	26-6101-50	500 pmols
Custom SPCT (Synthetic Positive Control Template); 131-180mer; 2 nmols	26-6101-02	2 nmols