



Product Specification Summary

Omni-Clean™ Gel DNA Purification Systems

| | |
|----------------|--|
| Catalog Number | 40-4111-75 |
| Product Name | OCC-1; 75 mL; Omni-Clean DNA Purification System |
| Size | 75 mL |
| Description | OCC-1; 75 mL; Omni-Clean DNA Purification System |
| Component/Note | OCC-1; 75 mL; Chaotrope blend 1 component of Omni-Clean DNA Purification System (Catalog#: 40-4110-XX) OCC-1; 75 mL; Chaotrope blend 1 component of Omni-Clean DNA Purification System (Catalog#: 40-4110-XX) |

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17356>



Product Manual <http://www.genelink.com/Literature/ps/M40-41XX-XX V3.2.pdf>



Related Products

| Product | Catalog No | Size |
|--|------------|-------------------|
| Omni-Clean™ Gel DNA Beads Purification System; 100 purifications | 40-4110-10 | 100 Purifications |
| Omni-Clean™ Gel DNA Beads Purification System; 500 purifications | 40-4110-50 | 500 Purifications |
| Omni-Clean™ Gel DNA Spin Column Purification System; 100 purifications | 40-4120-10 | 100 Purifications |
| Omni-Clean™ Gel DNA Spin Purification; 500 purifications | 40-4120-50 | 500 Purifications |
| OCC-1; 75 mL; Omni-Clean DNA Purification System | 40-4111-75 | 75 mL |
| Omni-Mag™ Lysis Buffer; 75 mL | 40-4119-75 | 75 mL |