

# Product Specification Summary

## Fluorescent Molecular Probes

Catalog Number 26-8124-02  
Product Name Fluorescent Probe, 5'-Cy 5/3'-Quencher; 15 nmol purified yield  
Size 15 nmols  
Description Fluorophore-Quencher Probe, 5'-Cy 5/3'-Quencher; 15 nmol purified yield  
Component/Note Fluorescent Probe, 5'-Cy 5/3'-Quencher; 15 nmol purified yield  
Scan the QR Code or visit the following links  
Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17423>

Product Manual [http://www.genelink.com/Literature/ps/Fluorescent\\_oligo\\_recons\\_V3.2.pdf](http://www.genelink.com/Literature/ps/Fluorescent_oligo_recons_V3.2.pdf)

Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>

### Related Products

Product	Catalog No	Size
SmartBase™ Fluorescent Molecular Probes	SB-FMP	1 each
QPCR_Manual	M_QPCR	1 each
Fluorescent Probe, 5'-Fam/3'-Quencher, 600 nmol purified yield	26-8121-10	600 nmols