



# Product Specification Summary

## GeneAssays™

Catalog Number	41-1102-01
Product Name	GeneAssays Custom qPCR Probe (NHS) & Primer mix, 1000 x 25uL reactions
Size	1000 Reactions
Description	GeneAssays Custom qPCR Probe (NHS) & Primer mix, 1000 x 25uL reactions
Component/Note	Custom probe with fluorophore that are NHS conjugated & primer mix as per customer supplied sequence.

Custom probe & primer mix as per customer supplied sequence and specifications. Price listed for AF Dyes. Prices vary for other dyes.  
GeneAssays Custom qPCR Probe & Primer mix, 1000 x 25uL reactions. Supplied dried in one tube as per customer specified primer and probe concentration.  
See 41-1101-01 for fluorophores that are directly coupled.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17426>



Product Manual [http://www.genelink.com/Literature/ps/GeneAssay\\_Custom\\_qPCR\\_Assay\\_PS\\_V3.2.pdf](http://www.genelink.com/Literature/ps/GeneAssay_Custom_qPCR_Assay_PS_V3.2.pdf)



Product MSDS [http://www.genelink.com/Literature/ps/MSDS\\_nonhaz\\_GeneAssays.pdf](http://www.genelink.com/Literature/ps/MSDS_nonhaz_GeneAssays.pdf)



### Related Products

Product	Catalog No	Size
Custom SPCT (Synthetic Positive Control Template); up to 130mer; 500 pmols	26-6100-50	500 pmol
Custom SPCT (Synthetic Positive Control Template); 131-180mer; 500 pmols	26-6101-50	500 pmols
Custom SPCT (Synthetic Positive Control Template); 131-180mer; 2 nmols	26-6101-02	2 nmols