



## Product Specification Summary

### Genemer™

Catalog Number	40-2036-10
Product Name	Factor VIII Genemer Pack (Hemophilia A); 10 nmols
Size	10 nmols
Description	Gene Detection System
Component/Note	Contains 27 Genemers
Component/Note	Genemer Pack

Genemer™ are optimized primer pairs for amplification of the gene fragment of interest, particularly those spanning a mutation. A wide range is available from single base pair mutation sites, the complex triple repeat disorders and various pathogens. Genemer™ kits are complete, easy-to-use systems for reliable genotyping of large repeats in certain triple repeat disorder. The Genemer™ product contains 10 nmols of the primer pair that is sufficient for 400 regular 50 µl amplifications. The Genemer™ Kit contains reagents as well as control DNA sufficient for 100 amplifications.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=207>



Product Manual

[http://www.genelink.com/Literature/ps/M40-2036-10\\_Ver3.1.pdf](http://www.genelink.com/Literature/ps/M40-2036-10_Ver3.1.pdf)



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
Cystic Fibrosis 7 Mutation Panel Genotyping Kit; 1 Kit	40-2029-07K	1 Kit
Fragile X Genemer Kit; 1 Kit	40-2004-11	1 Kit
DRPLA Genemer; 10 nmols	40-2042-10	10 nmols
HGH Genemer (Human Growth Hormone); 10 nmols	40-2024-10	10 nmols
STS Genemer (Steroid Sulfatase); 10 nmols	40-2023-10	10 nmols
Mouse-coxI Genemer; 10 nmols	40-1502-10	10 nmols
Rat-coxI Genemer; 10 nmols	40-1503-10	10 nmols
C-Hamster-Cyt.B Genemer; 10 nmols	40-1504-10	10 nmols
AG Monkey-coxI Genemer; 10 nmols	40-1505-10	10 nmols
Dog-coxI Genemer; 10 nmols	40-1506-10	10 nmols
Cat-coxI-For1; 10 nmols	40-1507-10F	10 nmols
Horse-coxI Genemer; 10 nmols	40-1509-10	10 nmols
Rabbit-coxI Genemer; 10 nmols	40-1508-10	10 nmols
Human-coxI Genemer; 10 nmols	40-1501-10	10 nmols
X alphoid repeat Genemer; 10 nmoles	40-2021-10	10 nmols