



# Product Specification Summary

## Genetic Tools & Reagents

Catalog Number	40-5025-10
Product Name	Maleic acid buffer 10X (Buffer M 10X); 100 mL
Size	100 mL
Description	Maleic acid buffer 10X (Buffer M 10X)
Component/Note	Maleic acid buffer 10X (Buffer M 10X)

Maleic acid buffer 10X (Buffer M 10X. Maleic acid buffer solution is used in non-radioactive digoxigenin based Southern blot analysis after stringent washing post hybridization with SSC/SDS based buffer solutions. Prior to anti-DIG-AP conjugate incubation the membrane is equilibrated in 1X Maleic after hybridization and then incubated in anti-DIG-AP conjugate diluted in Buffer MB composed of Maleic Acid buffer and Blocking reagent.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=4>



Product Manual [http://www.genelink.com/Literature/ps/Dig\\_Southern\\_Reagents\\_V3.4.pdf](http://www.genelink.com/Literature/ps/Dig_Southern_Reagents_V3.4.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Omni-Array™ mRNA Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit