



# Product Specification Summary

## RAPD Decamer Sets

Catalog Number	40-0001-15
Product Name	GL RAPD Decamer 15-ACAD Set
Size	1 set
Description	GL RAPD Decamer 15-ACAD Set
Component/Note	Supplied Lyophilized
Component/Note	40 Decamers 10 nmols each

Arbitrary sequence decamer sets for rapid discerning of DNA polymorphism. There are 30 sets available each containing 40 primers with a 60-70% GC content. These primers are used in amplification reactions for determining Random Amplified Polymorphic DNA (RAPD) fragments.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=735>



Product Manual <http://www.genelink.com/Literature/ps/PS40-0001-15.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
GL RAPD Decamer Set 1- 30	40-0001-00	1 set
GL RAPD Decamer 10-ST Set	40-0001-10	1 set