



# Product Specification Summary

## Genemer™

Catalog Number 40-2024-10  
Product Name HGH Genemer (Human Growth Hormone); 10 nmols  
Size 10 nmols

Genemer™ are optimized primer pairs for amplification of the gene fragment of interest, particularly those spanning a mutation. A wide range is available from single base pair mutation sites, the complex triple repeat disorders and various pathogens. Genemer™ kits are complete, easy-to-use systems for reliable genotyping of large repeats in certain triple repeat disorder. The Genemer™ product contains 10 nmols of the primer pair that is sufficient for 400 regular 50 µl amplifications. The Genemer™ Kit contains reagents as well as control DNA sufficient for 100 amplifications.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=79>



Product Manual <http://www.genelink.com/Literature/ps/ps40-2024-10.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Cystic Fibrosis 7 Mutation Panel Genotyping Kit; 1 Kit	40-2029-07K	1 Kit
Fragile X Genemer Kit; 1 Kit	40-2004-11	1 Kit
DRPLA Genemer; 10 nmols	40-2042-10	10 nmols