



# Product Specification Summary

## Genemer™ Kits

Catalog Number 40-2001-11K  
Product Name Sickle Cell Genemer Kit; 1 Kit  
Size 1 Kit  
Description Sickle Cell Genemer Kit; 1 Kit

The Sickle Cell Genemer kit is a ready to use kit with all components required to genotype HbS mutation. Simply add sample DNA. The kit is based on HbS mutation spanning fragment amplification followed by Dde I digestion and electrophoretic resolution.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=879>



Product Manual [http://www.genelink.com/Literature/ps/M40-2001-11K\\_V4.2.pdf](http://www.genelink.com/Literature/ps/M40-2001-11K_V4.2.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
SCA1 Genemer Kit; 1 Kit	40-2037-11	1 Kit
SCA2 Genemer Kit; 1 Kit	40-2038-11	1 Kit
SCA3 Genemer Kit; 1 Kit	40-2039-11	1 Kit