



Product Specification Summary

Genemer™ Kits

| | |
|----------------|--|
| Catalog Number | 40-2036-11 |
| Product Name | Factor VIII (Hemophilia) Genemer™ Kit; 1 Kit |
| Size | 1 Kit |
| Description | Factor VIII (Hemophilia) Genemer™ Kit |

Genemer™ are optimized primer pairs for amplification of the gene fragment of interest, particularly those spanning a mutation. A wide range is available from single base pair mutation sites, the complex triple repeat disorders and various pathogens. Genemer™ kits are complete, easy-to-use systems for reliable genotyping of large repeats in certain triple repeat disorder. The Genemer™ product contains 10 nmols of the primer pair that is sufficient for 400 regular 50 µl amplifications. The Genemer™ Kit contains reagents as well as control DNA sufficient for 100 amplifications.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=889>



Product Manual http://www.genelink.com/Literature/ps/M40-2036-10_Ver3.1.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|-------------------------|------------|-------|
| SCA1 Genemer Kit; 1 Kit | 40-2037-11 | 1 Kit |
| SCA2 Genemer Kit; 1 Kit | 40-2038-11 | 1 Kit |
| SCA3 Genemer Kit; 1 Kit | 40-2039-11 | |