



Product Specification Summary

Genetic Tools & Reagents

Catalog Number 40-3006-01
Product Name Omni- Marker™ Low unlabeled; 100 uL
Size 100 uL

Omni-Marker™ Universal and Low are unlabeled DNA markers containing a blend of fragments ranging from 50 base pairs to 10 kb. The universal contains fragments of the following sizes; 10 kb, 8 kb, 6 kb, 4 kb, 3 kb, 2 kb, 1.55, 1.4 kb, 1 kb, 750 bp, 500 bp, 400 bp, all the bands and “low” version contains fragments from 50 bp to 2kb. The low Omni-Marker™ is ideal for routine PCR gels. These are ready to load containing BPB/XC dyes. A loading of 5 µl is sufficient per lane.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=903>



Product Manual <http://www.genelink.com/Literature/ps/Omni-Marker ver2.3.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|--|------------|--------|
| Omni-Array™ mRN Antisense strand amplification kit | 08-0021-02 | 1 Kit |
| Omni-Array™ Sense Strand mRNA Amplification Kit | 08-0011-02 | 1 Kit |
| Omni-Marker™ Universal unlabeled; 100 uL | 40-3005-01 | 100 uL |