



Product Specification Summary

Genetic Tools & Reagents

Catalog Number 40-3072-05
Product Name Omni-Marker™ GScan-2 Tet labeled 50 bp- 600 bp
Size 500 ul

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=913>



Product Manual <http://www.genelink.com/Literature/ps/PS40-3062-XX.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|--|------------|--------|
| Omni-Array™ mRN Antisense strand amplification kit | 08-0021-02 | 1 Kit |
| Omni-Array™ Sense Strand mRNA Amplification Kit | 08-0011-02 | 1 Kit |
| Omni-Marker™ Universal unlabeled; 100 uL | 40-3005-01 | 100 uL |