



# Product Specification Summary

## Genetic Tools & Reagents

Catalog Number 40-3081-01  
Product Name GScan Marker Hex Labeled 50 bp - 600 bp; 100 uL  
Size 100 uL

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=914>



Product Manual <http://www.genelink.com/Literature/ps/Omni-Marker ver2.3.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Omni-Array™ mRN Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL