

## **Product Specification Summary**

## Omni-Pure™ Plasmid DNA Purification Systems

Catalog Number 40-4020-05

Product Name Omni-Pure™ Plasmid Purification System; 500 purifications

Size 500 Purifications

Description Omni-Pure™ Plasmid DNA Purification System

Facile and rapid purification of DNA from varied sources can be performed using the Omni-Pure™ series of DNA, RNA and plasmid purification systems. Each purification system has been formulated, optimized and designed to yield the highest purity available with the starting sample volume specially geared towards the desired downstream application. Routine mini-preps of plasmid extraction are made even easier with consistent performance. Purification can be performed with a maximum of 3 ml of cells yielding up to 20 µg of purified DNA. The convenient spin column method can be scaled up by using multiple columns and processed in less than 30 minutes. The purified DNA is of high quality suitable for all molecular biology applications including direct use in fluorescent automated sequencing methods.



## Scan the QR Code or visit the following links

Product Infomation http://www.genelink.com/geneprodsite/product.asp?p=989



Product Manual http://www.genelink.com/Literature/ps/ps40-4020-05.pdf



Product MSDS http://www.genelink.com/Literature/ps/MSDSNH.pdf



## **Related Products**

Product	Catalog No	Size
Omni-Pure™ Plasmid Purification System; 100 purifications	40-4020-01	100 Purifications
Omni-Pure™ Plasmid Purification System; 500 purifications	40-4020-05	500 Purifications



8 Westchester Plaza, Suite 130, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | www.genelink.com | SS40-4020-05.pdf Print Date Version : October 24, 2025 Page 2 of 2