



## Product Specification Summary

### Myotonic Dystrophy

Catalog Number	40-2026-40
Product Name	Myotonic Dystrophy GeneProber GLDM4 unlabeled Probe; 500 ng
Size	500 ng
Description	Myotonic Dystrophy Probe
Component/Note	Myotonic Dystrophy GeneProber GLDM4 unlabeled Probe; 500 ng
Storage	Store at -20°C

**Myotonic Dystrophy GeneProber GLDM4 unlabeled probe**

Myotonic dystrophy CTG triple repeat spanning region unlabeled probe for radioactive or chemiluminescent labeling and Southern blot detection of Bam HI, PstI, BglI and EcoRI digested DNA. Suitable for random primer labeling.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=631>



Product Manual

[http://www.genelink.com/Literature/ps/M40-2026-40\\_V6.2.pdf](http://www.genelink.com/Literature/ps/M40-2026-40_V6.2.pdf)



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
Myotonic Dystrophy GeneProber GLDM2 (Pst I) unlabeled Probe; 500 ng	40-2026-39	500 ng
Myotonic Dystrophy GScan™ TET Kit	40-2026-15TT	1 Kit
Myotonic Dystrophy GScan™ HEX Kit	40-2026-15HX	1 Kit
Myotonic Dystrophy GScan™ Cy3 Kit	40-2026-15C3	1 Kit
Myotonic Dystrophy GScan™ Cy5 Kit	40-2026-15C5	1 Kit
Myotonic Dystrophy GeneProber GLDM3 unlabeled Probe; 500 ng	40-2026-38	500 ng
Myotonic Dystrophy GScan™ 6-FAM Kit; 1 Kit	40-2026-15FM	1 Kit
Myotonic Dystrophy GScan™ 6-FAM Kit (20rxns); 1 Kit	40-2026-15FMS	1 Kit
Myotonic Dystrophy Genemer Kit; 1 Kit	40-2026-11	1 Kit
DM CTG-6-MB-CY3; 25 ug	45-0210-06	25 ug
DM CTG-6-MB-CY5; 25 ug	45-0220-06	25 ug
DM CTG-6-MB-FAM; 25 ug	45-0230-06	25 ug
DM CTG-6-MB-FI; 25 ug	45-0240-06	25 ug
DM CTG-6-MB-HEX; 25 ug	45-0250-06	25 ug
DM CTG-6-MB-TET; 25 ug	45-0260-06	25 ug